

2

7

PubChem » Substance Summary



Rhodamine 123 - Substance Summary

(SID: 50063483)

A fluorescent probe with low toxicity which is a potent substrate for P-glycoprotein and the bacterial multidrug efflux transporter. It is used to assess mitochondrial bioenergetics in living cells and to measure the efflux activity of P-glycoprotein in both normal and malignant cells. (Leukemia 1997;11(7):1124-30)

☐ Table of Contents

- Drug and Chemical Info
- Data Source
- Synonyms
- Properties
- Descriptors
- Substance Info
- Exports



Drug and Chemical Info from MeSH: (Total:1) 2

Rhodamine 123

Pharmacological Action

Fluorescent Dyes

Classification

Heterocyclic Compounds

Heterocyclic Compounds, 3-Ring

Xanthenes


Rhodamines

Rhodamine 123

PubMed Choose by Subheadings:

administration	adverse effects	analogs and	
and dosage	auverse effects	derivatives	
	antagonists and inhibitors	blood	The second of the second
chemistry	diagnostic use	isolation and a purification	20 00 0
metabolism	pharmacokinetics	pharmacology	200
therapeutic use	toxicity	urine	

PubMed MeSH Keyword Summary 2




■ Links

H-Bond Donor

H-Bond Acceptor 5

∞ 2



Rhodamine 123 EINECS 263-687-8
EINECS 263-687-8
RH 123
LS-162564
3,6-Diamino-9-(2-(methoxycarbonyl)phenyl)xanthylium chloride
Xanthylium, 3,6-diamino-9-(2-(methoxycarbonyl)phenyl)-, chloride
2-(6-Amino-3-imino-3H-xanthen-9-yl)benzoic acid methyl ester monohydrochloride
Benzoic acid, 2-(6-amino-3-imino-3H-xanthen-9-yl)-, methyl ester, monohydrochloride


Properties Computed from Structure: 2

	······································
Molecular Weight	380.82428 [g/mol]
Molecular Formula	$C_{21}H_{17}CIN_2O_3$
H-Bond Donor	2
H-Bond Acceptor	5
Rotatable Bond Count	3
Exact Mass	380.09277
MonoIsotopic Mass	380.09277
Topological Polar Surface Area	87.1
Heavy Atom Count	27
Formal Charge	0
Complexity	670
Isotope Atom Count	0
Defined Atom StereoCenter Count	.0
Undefined Atom StereoCenter Count	0
Defined Bond StereoCenter Count	0
Undefined Bond StereoCenter Coun	t0
Covalently-Bonded Unit Count	2

On Descriptors Computed from Structure: 2

IUPAC Name: [6-amino-9-(2-methoxycarbonylphenyl) xanthen-3-ylidene]azanium chloride
Canonical SMILES: COC(=0)
C1=CC=CC=C1C2=C3C=CC(=[NH2+])
C=C3OC4=C2C=CC(=C4)N.[Cl-]
InChI: InChI=1/C21H16N2O3.ClH/c1-25-21(24)15-5-3-2-4-14(15)20-16-8-6-12(22)10-18(16)26-19-11-13 (23)7-9-17(19)20;/h2-

11,22H,23H2,1H3;1H/fC21H17N2O3.Cl/h22H;1h/q+1;-1 ①





7

•

7

2





Properties Computed from Structure: 2

Molecular Weight	479.01034 [g/mol]
Molecular Formula	C ₂₈ H ₃₁ CIN ₂ O ₃
H-Bond Donor	2
H-Bond Acceptor	5
Rotatable Bond Count	7
Exact Mass	478.202321
MonoIsotopic Mass	478.202321
Topological Polar Surface Area	61.5
Heavy Atom Count	.34
Formal Charge	0
Complexity	823
Isotope Atom Count	0
Defined Atom StereoCenter Count	0
Undefined Atom StereoCenter Count	0/
Defined Bond StereoCenter Count	0
Undefined Bond StereoCenter Coun	t0
Covalently-Bonded Unit Count	2

🕍 Descriptors Computed from Structure: 🗵

IUPAC Name: [9-(2-ethoxycarbonylphenyl)-6ethylamino-2,7-dimethylxanthen-3-ylidene]-

ethylazanium chloride

Canonical SMILES: CCNC1=C(C=C2C(=C1)OC3=CC(= [NH+]CC)C(=CC3=C2C4=CC=CC=C4C(=0)OCC)C)C.

[CI-]

InChI: InChI=1/C28H30N2O3.ClH/c1-6-29-23-15-25-21(13-17(23)4)27(19-11-9-10-12-20(19)28(31)32-8-3)22-14-18(5)24(30-7-2)16-26(22)33-25;/h9-

16,29H,6-8H2,1-5H3;1H/b30-

24+;/fC28H31N2O3.Cl/h30H;1h/q+1;-1/b 2



Substance Info: 2

SID: 49854251 ☑ 🗐

Deposit Date: 2008-07-09 Modify Date: 2008-07-09

CID: 13806 🗹 🗷

Create Date: 2005-06-24 Parent CID: 13807 🗵 Unique Components: 2 Links

Related Substances: 2 Same: 14 Links

Similar Substances: 332 Links 2







PubChem » Substance Summary

HIDC iodide - Substance Summary (SID: 24895800)

Table of Contents


- Data Source
- Synonyms
- Properties
- Descriptors
- Substance Info
- Exports





Properties Computed from Structure: 2

Molecular Weight	510.45291 [g/mol]
Molecular Formula	$C_{27}H_{31}IN_2$
H-Bond Donor	0
H-Bond Acceptor	2
Rotatable Bond Count	3
Exact Mass	510.153192
MonoIsotopic Mass	510.153192
Topological Polar Surface Area	6.3
Heavy Atom Count	30
Formal Charge	0
Complexity	689
Isotope Atom Count	0
Defined Atom StereoCenter Count	0
Undefined Atom StereoCenter Count	0
Defined Bond StereoCenter Count	3
Undefined Bond StereoCenter Count	0
Covalently-Bonded Unit Count	2



∃ Links

Ť

Chemical Structure Search 2



Descriptors Computed from Structure: 🛭

IUPAC Name: 1,3,3-trimethyl-2-[(1E,3E,5Z)-5-(1,3,3trimethylindol-2-ylidene)penta-1,3-dienyl]indol-1-ium iodide Canonical SMILES: CC1(C2=CC=CC=C2[N+] (=C1C=CC=C3C(C4=CC=CC=C4N3C)(C)C)C)C.[I-]Isomeric SMILES: CC1(C2=CC=CC=C2[N+](=C1 C=CC=C/3C(C4=CC=CC=C4N3C)(C)C)C.[I-]InChI: InChI=1/C27H31N2.HI/c1-26(2)20-14-10-12-16-22(20)28(5)24(26)18-8-7-9-19-25-27(3,4)21-15-11-13-17-23(21)29(25)6;/h7-19H,1-6H3;1H/q+1;/p-1/fC27H31N2.I/h;1h/qm;-1 ②

🕽 Substance Info: 🛭

SID: 24895800 ② □

Deposit Date: 2007-07-16 Modify Date: 2008-01-21

CID: 6504806 ☑ ■

Create Date: 2006-05-01 Parent CID: 5477641 2 Unique Components: 2 Links


Related Substances: 2

Same: 3 Links

Same, Connectivity: 8 Links

Similar Substances: 787 Links 2





PubChem Substance PubChem | Go | PubMed | Entrez | Structure | PubChem | Help

PubChem » Substance Summary

DOTCI - Substance Summary (SID: 24859444)

Table of Contents

- Data Source
- Synonyms
- Properties
- Descriptors
- Substance Info
- Exports



Depositor: Sigma-Aldrich External ID: 32448_FLUKA

Depositor-Supplied Synonyms: (Total: 3) 🛭

Sort: Weight
DOTCI

32448_FLUKA 3,3'-Diethyloxatricarbocyanine iodide

Properties Computed from Structure: 2

	······································
Molecular Weight	512.38267 [g/mol]
Molecular Formula	$C_{25}H_{25}IN_2O_2$
H-Bond Donor	0
H-Bond Acceptor	4
Rotatable Bond Count	6
Exact Mass	512.096071
MonoIsotopic Mass	512.096071
Topological Polar Surface Area	29.5
Heavy Atom Count	30
Formal Charge	0
Complexity	601
Isotope Atom Count	0
Defined Atom StereoCenter Count	0
Undefined Atom StereoCenter	0
Count	
Defined Bond StereoCenter Count	4
Undefined Bond StereoCenter Coun	
Covalently-Bonded Unit Count	2


Compound Displayed PubChem

Compound ID 5702695	a
Molecular Weight 512.38267 [g/mol]	2
Molecular Formula C ₂₅ H ₂₅ IN ₂ O ₂	2
H-Bond Donor 0	2
H-Bond Acceptor 4	2

∃ Links

Chemical Structure Search 2

Descriptors Computed from Structure: 2



Related Substances: 🗵

Same: 5 Links

Same, Connectivity: 6 Links

Similar Substances: 898 Links 2



Oseroff



PubChem Substance PubMed | Entrez | Structure | PubChem | Help

PubChem » Substance Summary

DOTCI - Substance Summary (SID: 24859444)

Table of Contents

- Data Source
- Synonyms
- Properties
- Descriptors
- Substance Info
- Exports



Data Source: 🖼 🛭

Depositor: Sigma-Aldrich External ID: 32448_FLUKA


Depositor-Supplied Synonyms: (Total: 3) 🛭

Sort: Weight

DOTCI

32448 FLUKA

3,3'-Diethyloxatricarbocyanine iodide



Compound ID 5702695	(2)
Molecular Weight 512.38267 [g/mol]	2
Molecular Formula C ₂₅ H ₂₅ IN ₂ O ₂	Ø
H-Bond Donor 0	Ø
H-Bond Acceptor 4	2

Molecular Weight	512.38267 [g/mol]
Molecular Formula	$C_{25}H_{25}IN_2O_2$
H-Bond Donor	0
H-Bond Acceptor	4
Rotatable Bond Count	6
Exact Mass	512.096071
MonoIsotopic Mass	512.096071
Topological Polar Surface Area	29.5
Heavy Atom Count	30
Formal Charge	0
Complexity	601
Isotope Atom Count	0
Defined Atom StereoCenter Count	0
Undefined Atom StereoCenter Count	0
Defined Bond StereoCenter Count	4
Undefined Bond StereoCenter Count	0
Covalently-Bonded Unit Count	2

■ Links

Chemical Structure Search 2

Descriptors Computed from Structure: 2

Oseroff



PubChem Substance PubMed | Entrez | Structure | PubChem | Help

PubChem » Substance Summary



DODCI - Substance Summary

(SID: 661543 Version: 3) 2

RN given refers to iodide; RN for parent cpd not in Chemline 6/84

■ Table of Contents

- Drug and Chemical Info
- Data Source
- Synonyms
- Substance Info
- Exports



Drug and Chemical Info from MeSH: (Total:1) 2

3,3'-diethyloxadicarbocyanine



PubMed MeSH Keyword Summary 2

Data Source: 🗾 🛭

Depositor: ChemIDplus External ID: 014806509

Depositor-Supplied Synonyms: (Total: 12)

Display: Next 2 | All | Sort: Weight

DODCI

Eastman 11219

3,3'-Diethyloxadicarbocyanine • Diethyloxadicarbocyanine iodide

EINECS 238-873-7

3,3'-Diethyloxodicarbocyanine iodide

NSC 270081

3-Ethyl-2-(5-(3-ethyl-2(3H)-benzoxazolylidene)-1,3pentadienyl)benzoxazolium iodide

3-Ethyl-2-(5-(3-ethyl-3H-benzoxazol-2-ylidene)penta-


1,3-dienyl)benzoxazolium iodide

Benzoxazolium, 3-ethyl-2-(5-(3-ethyl-2(3H)-

benzoxazolylidene)-1,3-pentadienyl)-, iodide

Substance Info: 2





■ Links

NLM Toxicology (1) 2

Chemical Structure Search 2

SID: 661543 ② □ Deposit Date: 2005-08-08 Modify Date: 2006-04-29

ASN1 Display Display SDF Display SDF Save

Cyanine

You can support Wikipedia by making a tax-deductible donation.

From Wikipedia, the free encyclopedia

Cyanine is a non-systematic name of a synthetic dye family belonging to polymethine group.

Contents

- 1 Structure
- 2 Uses
- 3 Cy3 and Cy5
 - 3.1 Nomenclature and Structure
 - 3.2 Spectral characteristics
 - 3.3 Cye dye alternatives
 - 3.4 Cy5 ozone susceptibility
 - 3.5 Applications
 - 3.5.1 Nucleic acid labeling
 - 3.5.2 Protein labeling
- 4 References
- 5 See also

CH₃CH₂-NCH₃

Structure

There are three different types of cyanines:

Streptocyanines or open chain cyanines:

$$R_2N^+=CH[CH=CH]_n-NR_2$$
 (I)

Hemicyanines:


$$Aryl=N^{+}=CH[CH=CH]_{n}-NR_{2}$$
 (II)

Closed chain cyanines:^[1]

$$Aryl=N^+=CH[CH=CH]_n-N=Aryl$$
 (III)

where two nitrogens are joined by a polymethine chain^[2]. Both nitrogens are each independently part of a heteroaromatic moiety, such as pyrrole, imidazole, thiazole, pyridine, quinoline, indole, benzothiazole, etc.

Cyanines were first synthesized over a century ago, and there are a large number reported in the literature.



Uses

They have many uses as fluorescent dyes. Depending on the structure, they cover the spectrum from IR to UV.

They were originally used, and still are, to increase the sensitivity range of photographic emulsions, i.e., to increase the range of wavelengths which will form an image on the film.

They are used in CD-R and DVD-R media. The ones used are mostly green or light blue in color, and are chemically unstable. This makes cyanine discs unsuitable for archival CD and DVD use, as they can fade and become unreadable in a few years.

Cy3 and Cy5

Cy3 and Cy5 are reactive water-soluble fluorescent dyes of the cyanine dye family. Cy3 possesses red fluorescence, and Cy5 -

far-red. They are usually synthesized with reactive groups on either one or both of the nitrogen side chains so that they can be chemically linked to either nucleic acids or protein molecules. Labeling is done for visualization and quantification purposes. They are used in a wide variety of biological applications including comparative genomic hybridization and in gene chips, which are used in transcriptomics. They are also used to label proteins for various studies including proteomics.

Nomenclature and Structure

Standard chemical names specify exactly the chemical structure of the molecule. The Cy3 and Cy5 nomenclature was first proposed by Ernst, et al. [2] in 1989, and is non-standard, since it gives no hint of their chemical structures. In the original paper the number designated the count of the methines (as shown), and the side chains were unspecified. Thus various structures are designated Cy3 and Cy5 in the literature.

The R groups do not have to be identical. In the dyes as used they are short aliphatic chains one or both of which ends in a highly reactive moieties such as N-hydroxysuccinimide or maleimide.

Cy3

Spectral characteristics

Cy3 is excited maximally at 550 nm and emits maximally at 570 nm, in the red part of the spectrum; quantum yield is 0.15; FW=766.

Cy5 is excited maximally at 649 nm and emits maximally at 670 nm, in the far red part of the spectrum; quantum yield is 0.28. FW=792.

The scanners actually use different laser emission wavelengths (typically 532 nm and 635 nm) and filter wavelengths (550-600 nm and 655-695 nm) to avoid background contamination. They are thus able to easily distinguish between two samples when one sample has been labeled with Cy3 and the other labeled with Cy5. They are also able to quantitate the amount of labeling in either sample.

Cye dye alternatives

Alexa dyes and Dylight dyes can be used interchangeably with Cye dyes in most biochemical applications.

Cy5 ozone susceptibility

In 2003, researchers at Inpharmatics and Agilent reported in *Analytical Chemistry* that microarrays which used Cy5 were susceptible to intermittent data quality decrease caused by environmental ozone. Exposures to ozone levels above 5-10 ppb for 10-30 seconds were reported to decrease the reproducibility of Cy5 microarrays. Much higher levels of ozone (>100 ppb) were required to observe an effect in Cy3 [3].

Applications

Nucleic acid labeling

In microarray experiments DNA or RNA is labeled with either Cy3 or Cy5 that has been synthesized to carry an N-hydroxysuccinimidyl ester (NHS-ester) reactive group. Since, NHS-esters react readily only with aliphatic amine groups, which nucleic acids lack, nucleotides have to be modified with aminoallyl groups. This is done through incorporating aminoallyl-modified nucleotides during synthesis reactions. A good ratio is a label every 60 bases such that the labels are not too close to each other, thus resulting in quenching effects.

Protein labeling

For protein labeling, Cy3 and Cy5 dyes sometimes bear maleimide reactive groups instead. The maleimide functionality allows conjugation of the fluorescent dye to the sulfhydryl group of cysteine residues. Cysteines can be added and removed from the protein domain of interest via PCR mutagenesis.

Cy5, is sensitive to the electronic environment it resides in. Changes in the conformation of the protein it is attached to will

produce an enhancement or quenching of the emission. The rate of this change can be measured to determine enzyme kinetic parameters. The dyes can be used for similar purposes in FRET experiments.

Cy3 and Cy5 are used in proteomics experiments so that samples from two sources can be mixed and run together thorough the separation process^[4]. This eliminates variations due to differing experimental conditions that are inevitable if the samples were run separately. These variations make it extremely difficult, if not impossible, to use computers to automate the acquisition of the data after the separation is complete. Using these dyes makes the automation trivial.

References

- 1. ^ Johannes, H.H.: Cyanine: Direkte Funktionalisierung, Oligomerisierung, linear und nichtlinear optische Eigenschaften, Dissertation TU Braunschweig, 2000
- 2. ^ a b Ernst LA, Gupta RK, Mujumdar RB, Waggoner AS. Cyanine dye labeling reagents for sulfhydryl groups. Cytometry. 1989;10(1):3-10. PMID: 2917472
- 3. ^ Fare TL, Coffey EM, Hongyue D, et.al. Effects of Atmospheric Ozone on Microarray Data Quality. Analytical Chemistry. 2003;75:4672-4675. [1]
- 4. ^ Unlu M, Morgan ME, Minden JS. Difference gel electrophoresis: a single gel method for detecting changes in protein extracts. Electrophoresis. 1997;18(11):2071-7. PMID: 9420172

See also

- N-hydroxysuccinimide
- Maleimide

Retrieved from "http://en.wikipedia.org/wiki/Cyanine"
Categories: Dyes | Fluorescent dyes | Cyanine dyes | Quaternary ammonium compounds

- This page was last modified on 12 July 2008, at 08:36.
- All text is available under the terms of the GNU Free Documentation License. (See Copyrights for details.) Wikipedia® is a registered trademark of the Wikimedia Foundation, Inc., a U.S. registered 501(c)(3) tax-deductible nonprofit charity.